



Inverness[®]
SAFE EAR PIERCING

steps to ear piercing

step 1 Cleanse



Put on disposable gloves. Open cassette and remove alcohol prep pad. Using the alcohol prep pad, cleanse the front & back of the customer's ears. While ears are drying, provide customer with verbal aftercare instructions.


step 2 Mark



Mark ears with Inverness surgical marking pen. All marks should be placed on both ears prior to piercing. Always have customer approved placement of marks prior to piercing.

step 3

Load



Peel back each sterile cassette earring packing; Insert the earring cartridge into the instrument. You should hear it “snap” in.

To Aim & Pierce

Hold the Instrument firmly in one hand with the handle fitting in the palm of your hand. While keeping the back of the cartridge against the back of the ear, hold the ear with the other hand and place a finger against the back of the cartridge for support.

Apply firm/even pressure along the handle of the instrument and slowly depress the handle to push the inner capsule slightly out of the cartridge and towards the mark.

Aim carefully, making sure the capsule is lined up to mark & the instrument is level to the ear.

ATTENTION!!! DO NOT aim at the dot with the tip of the earring post! Forcing the earring post out of the capsule prior to piercing for the purpose of using the point of the earring to line up to the dot is called “TIPPING”, this action can result in misfires! Aim at the dot by using the end of the extended inner capsule only!

step

4

Pierce



Quickly & Firmly
depressing handle all of
the way until it can go no
further.

Hold . . . 1 . . . 2 . . . 3

Hold the handle to a count of 3 to ensure the plastic finger releases the earring back clutch

Move the instrument down and away carefully to make sure the piercing has been completed.


Repeat the above technique, pierce one ear and then the next.

Show the customer her/his newly pierced ears!

step

5

Eject



Eject the
cartridge and
dispose of the
trash.